

COLOUR CODING FOR IDENTIFICATION

Hypodermic Needle for Single Use - Colour Coding for Identification (As per ISO 6009:2016)		
Gauge (Size)	Nominal Outer Diameter Of Needle (in mm)	Colour Code
34	0.18	Orange
33	0.2	Black
32	0.23	Deep Green
31	0.25	White
30	0.3	Yellow
29	0.33	Red
28	0.36	Blue-Green
27	0.4	Medium Gray
26	0.45	Brown
25	0.5	Orange
24	0.55	Medium Purple
23	0.6	Deep Blue
22	0.7	Black
21	0.8	Deep Green
20	0.9	Yellow
19	1	Cream
18	1.2	Pink
17	1.4	Red-Violet
16	1.6	White
15	1.8	Blue-Gray
14	2.1	Pale Green
13	2.4	Purple
12	2.7	Pale Blue
11	3	Yellow-Green
10	3.4	Olive-Brown